

Surgical tie-on face masks

Irema has been manufacturing face masks under the Facemate® brand in Ireland for over 30 years.

All Facemate surgical masks are compliant with European Medical Device Regulation 2017/745 and tested by internationally accredited independent laboratories. They are manufactured using Irema's proprietary melt-blown filter material and meet the highest industry standards.

Features

- Latex and glass fibre free
- Lightweight
- 100% synthetic
- 3-ply or 4-ply
- Hypoallergenic materials do not cause skin irritation
- Lowest possible breathing resistance
- Available in blue and green

Specification

- EU Medical Device Regulation: 2017/745
- Test Standard: CE EN 14683
- CE EN149:2001 A1 2009
- Manufactured under ISO9001:2000 quality system and GMP CE marking approval
- Bacterial Filtration Efficiency BFE $\geq 98\%$
- Breathing resistance (mbar) 8.0 L/min: $<40.0 \text{ Pa/cm}^2$
- Type IIR splash resistant $>120 \text{ mmHG}$



Facemate surgical masks are packaged 50 masks per box (unless otherwise stated) with 6 boxes per case.

^ **Contact Sean Cooke**
scooke@irema.ie

Kilmallock Business Park
Kilmallock
Co Limerick
Ireland

email: info@irema.ie
Tel: + 353 63 98544

Facemate Anti-splash

Facemate Anti-splash Type IIR mask contains an additional splash resistant layer to protect the wearer against fluid splashes.

- 4-ply, Type IIR
- 40 per box
- Blue (green and white available on request)



Facemate Anti-fog

The Facemate Anti-fog is ideal for people who wear glasses or use optical equipment. It has a highly absorbent top binding layer that prevents glasses or surgical optics from fogging.

- Type II
- Highly absorbent top layer
- Blue (available in green and white on request)



Facemate Extra-large

Facemate Extra-large is a larger fitting mask with a wider open area to provide increased facial coverage. The mask has a foam top binding to provide a comfortable seal.

- Type II
- Suited for people with beards
- Blue (available in green and white on request)



Facemate XtraComfort

Facemate XtraComfort surgical mask has a special inner layer of non-woven material that feels extra smooth and perfect for sensitive skin. It keeps the face cool during extended periods of use.

- Also available in earloop
- Type II and IIR
- Blue, green, or white



Kilmallock Business Park
Kilmallock
Co Limerick
Ireland

email: info@irema.ie
Tel: + 353 63 98544